


Address

Jl. Utan Kayu Raya No. 104 - 105A
Jakarta 13120, INDONESIA

Contact

Phone : +62 21 8511 303
Fax : +62 21 8507 633, 8194342
Email : info@kharisma-utama.id
kharisma-utama@kharisma-utama.com
Web : www.kharisma-utama.id



FURNITURE CATALOGUE





MODERN HOSPITAL FURNITURE

New level of comfort and safety

WITH SUPER QUALITY AND AFFORDABLE PRICE

MADE IN INDONESIA
KHARISMA UTAMA

“OUR HOSPITAL FURNITURE,
IS NOT JUST ABOUT FURNITURE - TO ENSURE AN
APPROPRIATE AND FUNCTIONAL WORKING
ENVIRONMENT FOR MEDICAL STAFF AND PATIENT,
IT MUST BE VERY COMFORTABLE AND WELL
INTEGRATED.




GUARANTEES HIGH QUALITY, EASY
AND FAST DISINFECTION AND
ERGONOMICS.

”

**Hospital furniture
Category**

004-01HB MCA Hospital Electric Bed *Superior 4*



Rp. 76.000.000,-



Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (central brake system)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

LKPP



SCAN ME




Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

004-01HB MCA Hospital Electric Bed *Superior 3*



Rp. 68.000.000,-



Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (central brake system)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

LKPP



SCAN ME




Certified Company
Cert No :44 221 127768




Certified Company
Cert No :16 00 B 14030

004-01HB MCA Hospital Electric Bed *Superior 2*



Rp. 59.000.000,-



Equipped with a **remote control** to adjust the bed position



The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



The **Base** made of strong and lightweight material



Side guards for patient safety are designed to be easy to operate

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

LKPP



SCAN ME




Certified Company
Cert No :44 221 127768




Certified Company
Cert No :16 00 B 14030

004-01HB MCA Hospital Electric Bed *Superior 1*



Rp. 49.124.000,-



Equipped with a **remote control** to adjust the bed position



The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



The **Base** made of strong and lightweight material



Side guards for patient safety are designed to be easy to operate

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

LKPP



SCAN ME




Certified Company
Cert No :44 221 127768




Certified Company
Cert No :16 00 B 14030

004-01HB

MCA Hospital Electric Bed with Latex Mattress

Superior 4

Rp. 82.500.000,-

Latex


Spesification :


- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (central brake system)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP

SCAN ME


004-01HB

MCA Hospital Electric Bed with Latex Mattress

Superior 3

Rp. 73.800.000,-

Latex


Spesification :


- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (central brake system)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP


SCAN ME


004-01HB

MCA Hospital Electric Bed with Latex Mattress

Superior 2


Rp. 65.000.000,-


Latex


Equipped with a **remote control** to adjust the bed position



The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



The **Base** made of strong and lightweight material



Side guards for patient safety are designed to be easy to operate

Spesification :

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP


SCAN ME




Certified Company
Cert No :44 221 127768





Certified Company
Cert No :16 00 B 14030

004-01HB

MCA Hospital Electric Bed with Latex Mattress

Superior 1


Rp. 55.400.000,-


Latex


Equipped with a **remote control** to adjust the bed position



The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



The **Base** made of strong and lightweight material



Side guards for patient safety are designed to be easy to operate

Spesification :

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP


SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

001-01HB

MCA Hospital Bed



Rp. 17.911.000,-



equipped with a **folding handle** to adjust the position of the bed, crank handle made of plastic or aluminum

The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



Side guards for patient safety are designed to be easy to operate



*note : base on stock in factory

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 60 cm

LKPP



SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

002-01HB

MCA Hospital Bed



Rp. 20.391.000,-



equipped with a **folding handle** to adjust the position of the bed, crank handle made of plastic or aluminum

The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



Side guards for patient safety are designed to be easy to operate



*note : base on stock in factory

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 60 cm

LKPP



SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

001-01HB

MCA Hospital Bed with Mattress



Rp. 20.300.000,-



*note : base on stock in factory

Spesification :

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 60 cm
- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP



SCAN ME

Certified Company
Cert No : 44 221 127768Certified Company
Cert No : 16 00 B 14030

002-01HB

MCA Hospital Bed with Mattress



Rp. 22.700.000,-



*note : base on stock in factory

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 60 cm
- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP



SCAN ME

Certified Company
Cert No : 44 221 127768Certified Company
Cert No : 16 00 B 14030

003-01HB

MCA Hospital Bed *Superior*



Rp. 28.290.000,-



equipped with a **folding handle** to adjust the position of the bed, crank handle made of plastic or aluminum

*note : base on stock in factory

The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



Side guards for patient safety are designed to be easy to operate



The **Base** made of strong and lightweight material

Spesification :

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Hi-low can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

LKPP


SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

011-01CB

MCA Children Bed



Rp. 20.207.000,-

LKPP


SCAN ME



Spesification :

- Main component made from steel and stainless steel pipe
- Back rise can be adjust manually
- Side guard can be set manually
- Equipped with 4 casters (2 caster with brake)
- Dimension :
± (172 Length x width 80) cm



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

MCA Children Bed with Mattres



SCAN ME

LKPP

Rp. 22.200.000,-



Spesification :


- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm



Foam

004-01HB


MCA Hospital Electric Bed with Super Foam Mattress

Superior 1


Rp. 52.028.000,-



Foam



Equipped with a **remote control** to adjust the bed position



The **head and foot boards** are made of ABS which is strong and easy to clean, easy to install and remove.



The **Base** made of strong and lightweight material



Side guards for patient safety are designed to be easy to operate

Spesification :

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP



002-01HB

MCA Hospital Bed



Rp. 39.500.000,-



latex

SPESIFICATION

002-01HB 2 Crank Hospital Bed

- Construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 60 cm

MCA 002-MTS Latex Foam Mattress

- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

MCA 003-01IS Infusion Stand

- Feature : Adjustable height with locking screw
- Upper stand : Stainless steel pipe and rod
- Lower stand : Stainless steel pipe
- Foot : From plastic star with 5 casters
- Dimension : ± Height 120 - 220 cm

MCA 003-05BC Bedside Cabinet

- Construction made from steel sheet
- Top table with guard railing made from stainless steel
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 1 rack shoes
- Equipped with 4 casters
- Dimension :
± Length 45 x Width 40 x Height 85 cm

MCA 001-03OBT Overbed Table

- Top table made from plywood laminated with formica
- Construction made from steel
- Height stand can be adjusted up and down easily use knob star
- Equipped with 4 casters
- Dimensions :
± Length 70 x Width 40 x Height 80 - 120 cm

LKPP



003-01HB

MCA Hospital Bed

Deluxe

Rp. 30.400.000,-

Side guards for patient safety are designed to be easy to operate

OPSI 1

OPSI 2

OPSI 3


The head and foot boards are made of ABS which is strong and easy to clean, easy to install and remove.

The Base made of strong and lightweight material

equipped with a **folding handle** to adjust the position of the bed, crank handle made of plastic or aluminum

*note : base on stock in factory

Spesification :

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Hi-low can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

LKPP

SCAN ME


Certified Company
Cert No :44 221 127768




Certified Company
Cert No :16 00 B 14030

004-01IB

MCA ICU Bed with Latex Mattress

Superior 1

Rp. 96.500.000,-

Latex


Equipped with a **remote control** to adjust the bed position

The head and foot boards are made of ABS which is strong and easy to clean, easy to install and remove.



The Base made of strong and lightweight material

Side guards for patient safety are designed to be easy to operate



Wheels with central brake

Spesification :

- Constructions made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by electric system use remote
- Knee raise can be adjusted by electric system use remote
- Hi-low can be adjusted by electric system use remote
- Trendelenburg can be adjusted by electric system use remote
- Power supply is 220V - 240V, 50/60Hz
- Capacity is 144 Watt
- Product weight is 105 kg
- Maximum patient weight is 150 kg
- Safe Workload is 180 kg
- Equipped 4 casters
- Dimensions :
± (Length 215 x Width 98 x Height 60) cm


SCAN ME

003-01HB

MCA Hospital Bed

Deluxe II


Rp. 37.000.000,-

Side guards for patient safety are designed to be easy to operate

OPSI 1

OPSI 2

OPSI 3


SCAN ME

LKPP
Foam


The head and foot boards are made of ABS which is strong and easy to clean, easy to install and remove.



The Base made of strong and lightweight material

equipped with a **folding handle** to adjust the position of the bed, crank handle made of plastic or aluminum


*note : base on stock in factory

Spesification :

- Bed** construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Hi-low can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm
- Mattress** main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm
- Bedside cabinet** top table made from ABS
- Retractable made from ABS
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 retractable
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 4 casters
- Dimension :
± Length 46 width 44 x Height 77 cm

004-01HB

MCA Electric Hospital Bed

Grand Luxury


Rp. 96.000.000,-


Latex

SPESIFICATION

004-01HB Electric Hospital Bed

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted by remote control
- Knee raise can be adjusted by remote control
- Hi-low can be adjusted by remote control
- Equipped 4 casters (central brake system)
- Dimension :
± Length 215 x Width 98 x Height 50-80 cm

MCA 002-MTS Latex Foam Mattress

- Mattress main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm

MCA 003-01IS Infusion Stand

- Feature : Adjustable height with locking screw
- Upper stand : Stainless steel pipe and rod
- Lower stand : Stainless steel pipe
- Foot : From plastic star with 5 casters
- Dimension : ± Height 120 - 220 cm

MCA 003-10BC Bedside Cabinet

- Top table made from ABS
- Retractable made from ABS
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 retractable
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 4 casters
- Dimension :
± Length 46 width 44 x Height 77 cm

MCA 001-02OBT Overbed Table

- Top table made from plywood laminated with formica
- Construction made from stainless steel
- Height stand can be adjusted up and down easily use knob star
- Equipped with 4 casters
- Dimensions :
± Length 70 x Width 40 x Height 80 - 120 cm

LKPP


SCAN ME



003-01HB MCA Hospital Bed *Deluxe III*



Rp. 45.500.000,-

Side guards for patient safety are designed to be easy to operate

OPSI 1



OPSI 2



OPSI 3



SPESIFIKATION

003-01HB Hospital Bed

- Bed construction made from steel pipe
- Mattress base made from steel sheet
- Back raise can be adjusted manually by crank
- Knee raise can be adjusted manually by crank
- Hi-low can be adjusted manually by crank
- Equipped 4 casters (2 casters with brakes)
- Dimension :
- $\pm$ Length 215 x Width 98 x Height 50-80 cm

MCA 001-MTS Foam Mattress

- Mattress main material made from super foam and covered with synthetic leather
- Size of product :
- $\pm$ Length 200 x Width 90 x Height 10 cm

MCA 003-01IS Infusion Stand

- Feature : Adjustable height with locking screw
- Upper stand : Stainless steel pipe and rod
- Lower stand : Stainless steel pipe
- Foot : From plastic star with 5 casters
- Dimension : $\pm$ Height 120 - 220 cm

MCA 003-10BC Bedside Cabinet

- Top table made from ABS
- Retractable made from ABS
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 retractable
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 4 casters
- Dimension :
- $\pm$ Length 46 width 44 x Height 77 cm

MCA 001-02OBT Overbed Table

- Top table made from plywood laminated with formica
- Construction made from stainless steel
- Height stand can be adjusted up and down easily use knob star
- Equipped with 4 casters
- Dimensions :
- $\pm$ Length 70 x Width 40 x Height 80 - 120 cm

013-01GET, MCA Gynecological Examination Table



Rp.19.279.000,-

Spesification :

- Constructions from pipe
- Back raise can be set by manually system
- Leg section can be set by manually system
- Equipped with infusion stand, basin, knee support
- Dimensions :
- $\pm$ (Lenght 195 x Width 60) cm



SCAN ME

LKPP

021-01VB, MCA Verlos Bed

Rp. 21.318.000,-



Spesification :

- Constructions made from pipe
- Back raise can be set manually system
- Equipped with infusion stand, basin, knee support
- Equipped 4 casters when 2 casters with brakes
- Dimensions :
- $\pm$ (Lenght 210 x Width 80) cm



SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

001-01ES, MCA Emergency Strecher

Rp. 15.835.000,-



Spesification :

- Construction made from stainless steel pipe
- Infusion Stand pole made from stainless steel pipe
- Side Rail pole made from stainless steel pipe
- Mattress made from foam with vinyl cover
- Back Raise can be set by manually system
- Equipped with 4 casters (2 with brakes)
- Dimension :
± Length 200 x Width 60 x Height 85 cm

LKPP



SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

001-01HS, MCA Hydraulic Strecher

Rp. 69.990.000,-



Spesification :

- Construction made from stainless steel pipe, steel pipe, stainless steel sheet and steel sheet
- Head and foot panel, side guard, Oxygen rack
- Mattress made from foam with vinyl cover
- Back Raise adjusted by gas spring
- Hi-Low adjusted by hydraulic
- Equipped with 4 casters
- Dimension :
± Length 205 x Width 68 x Height 53 - 86 cm

Side Option



*note : base on stock in factory

LKPP



SCAN ME



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

011-01BB, MCA Baby Basket *Acrylic*



Spesification :

- Main component made from stainless steel pipe
- Basket made from acrylic
- Basket can be removable and tiltable
- Equipped with 4 casters
- Dimension :
 - Trolley
± (70 Length x 40 Width x 65 Height) cm
 - Basket
± (65 Length x 35 Width x 35 Height) cm

LKPP

SCAN ME

Rp. 7.991.000,-


011-01BB, MCA Baby Basket with Mattress *Acrylic*



Spesification :

- Bed main component made from stainless steel pipe
- Basket made from acrylic
- Basket can be removable and tiltable
- Equipped with 4 casters
- Dimension :
 - Trolley
± (70 Length x 40 Width x 65 Height) cm
 - Basket
± (65 Length x 35 Width x 35 Height) cm
- Main material made from super foam and covered with synthetic leather

LKPP

SCAN ME

Rp. 8.700.000,-


011-02BB, MCA Baby Basket *Plastic*



Spesification :

- Main component made from stainless steel pipe
- Basket made from plastic
- Basket can be removable and tiltable
- Equipped with 4 casters
- Dimension :
± (81 Length x 46 Width x 85 Height) cm

LKPP

SCAN ME

Rp. 8.267.000,-


011-02BB, MCA Baby Basket with Mattres *Plastic*



Spesification :

- Main component made from stainless steel pipe
- Basket made from plastic
- Basket can be removable and tiltable
- Equipped with 4 casters
- Dimension :
± (81 Length x 46 Width x 85 Height) cm
- Main material made from super foam and covered with synthetic leather

LKPP

SCAN ME

Rp. 8.950.000,-

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

003-01IS MCA Infusion Stand


LKPP


SCAN ME

Rp. 1.997.000,-

Spesification :

- Feature : Adjustable height with locking screw
- Upper stand : Stainless steel pipe and rod
- Lower stand : Stainless steel pipe
- Foot : From plastic star with 5 casters
- Dimension :
± Height 120 - 220 cm


Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

011-01ET MCA Examination Table


LKPP


SCAN ME

Rp. 9.277.000,-

Spesification :

- Construction made from stainless steel pipe
- Mattress made from foam with vinyl cover
- Equipped 4 legs with ruber foot
- Dimension :
± Length 200 x Width 60 x Height 80 cm

001-01RT, MCA Resuscitation Table


LKPP


SCAN ME

Rp. 15.065.000,-

Spesification :

- Main component made from steel
- Intended for baby
- Equipped with 2 drawer
- Equipped with cabinet locked system
- Dimension :
±Length 100 x Width 60 x Height 74 cm


Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

011-02ET MCA Examination Table


LKPP


SCAN ME

Rp. 14.400.000,-

Spesification :

- Construction made from stainless steel pipe
- Mattress made from foam with vinyl cover
- Back raise can be adjusted manually by gass spring
- Equipped 2 foot list
- Max patient weight is ± 120 kg
- Equipped 4 legs with ruber foot
- Dimension :
± Length 200 x Width 60 x Height 80 cm

MCA 002-04BC, Bedside Cabinet


LKPP

SCAN ME

Rp. 6.171.000,-

Spesification :

- Construction made from steel sheet
- Top table made from plywood laminated with formica
- Retractable made from plywood laminated with formica
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 drawer
- Equipped with 1 retractable
- Equipped with 1 cabinet
- Equipped with 1 rack shoes
- Equipped with 4 casters
- Dimension :

± Length 45 x Width 40 x Height 85 cm

MCA 002-06BC Bedside Cabinet


LKPP

SCAN ME

Rp. 7.165.000,-

Spesification :

- Construction made from steel sheet
- Top table with guard railing made from stainless steel
- Retractable made from plywood laminated with formica
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 drawer
- Equipped with 1 retractable
- Equipped with 1 cabinet
- Equipped with 1 rack shoes
- Equipped with 4 casters
- Dimensions :

± Length 45 x Width 40 x Height 85 cm

MCA 002-05BC, Bedside Cabinet


LKPP

SCAN ME

Rp. 5.695.000,-

Spesification :

- Construction made from steel sheet
- Top table with guard railing made from stainless steel
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 1 rack shoes
- Equipped with 4 casters
- Dimension :

± Length 45 x Width 40 x Height 85 cm

MCA 003-10BC Bedside Cabinet


LKPP

SCAN ME

Rp. 7.348.000,-

Spesification :

- Top table made from ABS
- Retractable made from ABS
- Towel hanger made from stainless steel
- Equipped with 1 towel hanger
- Equipped with 1 retractable
- Equipped with 1 drawer
- Equipped with 1 cabinet
- Equipped with 4 casters
- Dimension :

± Length 46 width 44 x Height 77 cm

MCA 002-05BS Bed Screen



LKPP

SCAN ME

Rp. 3.674.000,-

Spesification :

- Main material made from steel pipe
- Cloth can be detachable for cleaning
- 2 part attached together can be fold
- Equipped with 1 caster
- Dimension :
± Length (2 x 90) x Height 170 cm



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

MCA 002-06BS Bed Screen



LKPP

SCAN ME

Rp. 6.154.000,-

Spesification :

- Main material made from steel pipe
- Cloth can be detachable for cleaning
- 3 part attached together can be fold
- Equipped with 1 caster
- Dimension :
± Length (3 x 90) x Height 170 cm

MCA 001-04BST, Bowl Stand



LKPP

SCAN ME

Rp. 3.858.000,-

Spesification :

- Construction made from stainless steel pipe
- Bowl made from stainless steel
- Equipped with 2 bowl
- Equipped with 2 towel hanger
- Equipped with 5 casters (Star type)
- Dimension :
± Length 75 x Height 76 cm

MCA 001-01BDT, Baby Dressing Table



LKPP

SCAN ME

Rp. 15.065.000,-

Spesification :

- Main component made from steel
- Intended for baby
- Equipped with 2 drawer
- Equipped with cabinet locked system
- Dimension :
± Length 100 x Width 60 x Height 74 cm

MCA 001-01DT Dressing Trolley


LKPP

SCAN ME

Rp. 15.431.000,-

Spesification :

- Construction made from stainless steel pipe
- Equipped with 4 casters
- Equipped with 2 shelves
- Equipped with 1 handle
- Equipped with 1 waskom
- Equipped with 1 Bucket
- Dimension :
± Length 75x Width 50 x Height 85 cm


Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

MCA 001-02FS Foot Step


LKPP

SCAN ME

Rp. 1.906.000,-

Spesification :

- Main component made from stainless steel pipe
- Base made from plywood
- Have 2 step
- Dimension :
± Length 50 x Width 30 x Height 40 cm

MCA 001-01LHT, Linen Hamper Trolley


LKPP

SCAN ME

Rp. 46.017.000,-

Spesification :

- Construction made from stainless steel sheet and stainless steel pipe
- Equipped with guard railing
- Equipped with handle
- Equipped with 1 bracket hanger
- Equipped with 1 linen bag for dirty cloth
- Equipped with 1 cupboard with door lock system including 3 shelves for clean cloth
- Equipped with 4 casters (2 casters with brakes)
- Dimension :
± Length 108 x 63 Width x Height 107 cm

MCA 002-03LHT, Linen Hamper Trolley


LKPP

SCAN ME

Rp. 22.044.000,-

Spesification :

- Construction made from stainless steel sheet and stainless steel pipe
- Equipped with guard railing
- Equipped with handle
- Equipped with 1 bracket hanger
- Equipped with 1 linen bag for dirty cloth
- Equipped with 1 cupboard with door lock system including 3 shelves for clean cloth
- Equipped with 4 casters (2 casters with brakes)
- Dimension :
± Length 108 x 63 Width x Height 107 cm

MCA 001-01IC Instrument Cabinet



Spesification :

- Construction made from stainless steel sheet
- Equipped with 1 door with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 65 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 28.015.000,-

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

MCA 001-02IC Instrument Cabinet



Spesification :

- Construction made from stainless steel sheet
- Equipped with 2 door with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 120 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 37.108.000,-

LKPP

SCAN ME

Rp. 14.421.000,-

MCA 002-03IC Instrument Cabinet



Spesification :

- Construction made from steel sheet
- Equipped with 1 door with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 65 x Width 45 x Height 170 cm

MCA 002-04IC Instrument Cabinet



Spesification :

- Construction made from steel sheet
- Equipped with 2 door with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 120 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 15.431.000,-

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

MCA 001-01IT Instrument Trolley



Spesification :

- Construction made from stainless steel pipe and stainless sheet
- Equipped with 4 casters
- Equipped with 2 shelves
- Equipped with 1 handle
- Dimension :
± Length 60 x Width 40 x Height 80 cm

LKPP

SCAN ME

Rp. 5.162.000,-

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

MCA 001-04IT Instrument Trolley



Spesification :

- Construction made from stainless steel and stainless sheet
- Equipped with 4 casters
- Equipped with 3 shelves
- Dimension :
± Length 60 x Width 40 x Height 80 cm

LKPP

SCAN ME

Rp. 6.614.000,-

MCA 001-06IT Instrument Trolley



Spesification :

- Construction made from stainless steel and stainless sheet
- Drawer made from plastic
- Equipped with 4 casters
- Equipped with 2 shelves
- Equipped with 2 drawer
- Dimension :
± Length 60 x Width 40 x Height 80 cm

LKPP

SCAN ME

Rp. 5.934.000,-

Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

MCA 001-09IT Instrument Trolley



Spesification :

- Construction made from stainless steel and stainless sheet
- Equipped with 4 casters (2 with brake)
- Equipped with 2 shelves
- Equipped with 2 drawer with handle
- Dimension :
± Length 60 x Width 40 x Height 80 cm

LKPP

SCAN ME

Rp. 22.687.000,-

MCA 001-11IT Instrument Trolley



Spesification :

- Construction made from stainless steel and stainless sheet
- Drawer made from plastic
- Equipped with 4 casters (2 casters with brake)
- Equipped with 2 shelves
- Equipped with 2 drawer
- Dimension :
± Length 60 x Width 40 x Height 95 cm


LKPP

SCAN ME

Rp. 11.665.000,-

MCA 001-MTS Mattress Foam



Spesification :

- Main material made from super foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm


Certified Company
Cert No :44 221 127768

Certified Company
Cert No :16 00 B 14030

LKPP

SCAN ME

Rp. 2.904.000,-

MCA 002-03LC Laundry Cart



Spesification :

- Consturction made from steel
- Equipped with 4 caster (2 casters with brakes)
- Laundry bag made from linen cloth
- Equipped with 2 handle
- Dimension :
± Length 90 x Width 60 x Height 75 cm


LKPP

SCAN ME

Rp. 6.063.000,-

MCA 002-MTS Mattress Foam Latex



Spesification :

- Main material made from latex foam and covered with synthetic leather
- Size of product :
± Length 200 x Width 90 x Height 10 cm


LKPP

SCAN ME

Rp. 7.714.000,-

MCA 001-02MT, Mayo Table



LKPP

SCAN ME

Rp. 6.063.000,-

Spesification :

- Tray made from stainless steel sheet
- Main component made from stainless steel
- Height stand can be adjusted up and down easily use knob star
- Have 4 casters
- Dimensions :
± Length 60 x Width 40 x Height 80 - 120 cm

MCA 002-03OBT Overbed Table



LKPP


SCAN ME

Rp. 3.858.000,-

Spesification :

- Top table made from plywood laminated with formica
- Construction made from steel and stainless steel
- Height stand can be adjusted up and down easily use knob star
- Equipped with 4 casters
- Dimensions :
± Length 70 x Width 40 x Height 80 - 120 cm

MCA 001-02OBT Overbed Table



LKPP

SCAN ME

Rp. 5.916.000,-

Spesification :

- Top table made from plywood laminated with formica
- Construction made from stainless steel
- Height stand can be adjusted up and down easily use knob star
- Equipped with 4 casters
- Dimensions :
± Length 70 x Width 40 x Height 80 - 120 cm



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

MCA
Powered by **Kharisma Utama**

MCA 001-01MC Medicine Cabinet



Spesification :

- Construction made from stainless steel sheet
- Equipped with 1 door with lock and handle
- Equipped with 1 locker with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 65 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 24.947.000,-



Certified Company
Cert No :44.221.127768



Certified Company
Cert No :16.00.B.14030

MCA 001-02MC Medicine Cabinet



Spesification :

- Construction made from stainless steel sheet
- Equipped with 2 door with lock and handle
- Equipped with 2 locker with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 120 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 37.659.000,-

MCA 002-03MC Medicine Cabinet



Spesification :

- Construction made from steel sheet
- Equipped with 1 door with lock and handle
- Equipped with 1 locker with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 65 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 11.941.000,-

MCA 002-04MC Medicine Cabinet



Spesification :

- Construction made from steel sheet
- Equipped with 2 door with lock and handle
- Equipped with 2 locker with lock and handle
- Equipped with 3 shelves
- Dimension :
± Length 120 x Width 45 x Height 170 cm

LKPP

SCAN ME

Rp. 19.013.000,-

Terms and Conditions

Syarat dan Ketentuan



1. Pembelian dapat dilakukan melalui website www.lkpp.go.id
2. harga yang tercantum sudah termasuk PPn 11%.
3. Harga yang tercantum sudah termasuk 1x instalasi dan 1x training, apabila di perlukan training kembali maka biaya akan di kenakan oleh pembeli.
4. Harga yang tercatum sudah termasuk garansi selama 1 tahun, sesuai dengan ketentuan umum dan persyaratan garansi yang berlaku.
5. Harga yang tercantum belum termasuk ongkos kirim, yang akan ditentukan sesuai ongkos kirim actual yang berlaku saat pelaksanaan pengiriman di tempat dan lokasi customer.

