


TONSILLECTOMY SURGERY SET



HB 0682
Negus hemostatic forceps,
Fig.2, 19cm, 7 1/2"



HB 1852
Metzenbaum dissecting
scissors, curved, Fig.1,
18cm, 7"



HB 2680
Yankauer suction tube,
6mm, 30cm, 10 3/4"




HB 5711
Barnhill adenoid curette,
Fig.1, 21cm, 8 1/4"



HB 5712
Barnhill adenoid curette,
Fig.2, 21cm, 8 1/4"



HB 5713
Barnhill adenoid curette,
Fig.3, 21cm, 8 1/4"



HB 5830
Jennings mouth gag, 9cm, 3 1/2"



HB 5831
Jennings mouth gag, 11cm, 4 1/2"



HB 5891
Tobold tongue depressor,
solid, 14,5cm, 5 3/4"



HB 6004
Schnidt tonsil seizing forceps,
Fig.2, 19cm, 7 1/2"



HB 6005
White-Lillie tonsil seizing
forceps, 19cm, 7 1/2"



HB 6006
Blohmke tonsil seizing
forceps, 20cm, 8"

TONSILLECTOMY SURGERY SET



HB 5839
Davis-Boyle Set, complete with :
handle HB 5840-02, 2 silicon inserts
5 interchangeable tongue depressors HB 5841 – HB 5845



HB 6012
Eves tonsil snare,
28cm, 11"



HB 6018
Snare wires, cut,
(pack of 12 pcs.), 0,3x100mm



Fig. 1 Fig. 2 Fig. 3




HB 6027
Sluder-Ballenger tonsillectome
with tip Fig.1,2+3



HB 6019-03
Fowler tonsil knife,
20cm, 8"



HB 6021
Hurd tonsil dissector and
retractor, 23cm, 9"



HB 6022-01
Henke tonsil dissector,
Fig.1, 23cm, 9"